



Sanitary installation		MEIKO	In the building
B	Cold Water	1/2" NPT "Y" Strainer	Inlet pipe 3/4". Stop-cocks, 41°-77°F Flow: flow pressure >4.76 Gal/Minute., >14.5 PSI
C	Hot water	1/2" NPT "Y" Strainer	Inlet pipe 3/4". Stop-cocks, 110°-140°F Flow: flow pressure >4.76 Gal/Minute., >14.5 PSI
E	Waste Outlet TopLine40	3" "P" Trap	Waste pipe 3" DFU 4
F	Vent TopLine40	2" Vent Connection	Local Area Vent 2" 12.4 CFM
Electrical installation			Wall socket and main fuse. Protection: 3x15A. Protective measures as well as connection of the equi-potential bonding must be carried out according to IEC 364-5-54 and conform to the local and national power utilities regulations.
A	3 Ph ~ 208/230V, 60Hz, 4.2KW/4.6KW		
Protection class: Equipment Class 1, Pollution Degree 2, Installation Category II, Main switch,			

Connection dimensions are subject to tolerances.
It is neither allowed to give this drawing without our permissions to a third party for information or reproduction,
nor may it be used for competition purposes. All rights are reserved! Subject to technical changes.



MEIKO USA INC.
1349 Heil Quaker Blvd. LaVergne TN 37086
UNITED STATES OF AMERICA
TEL. (615)399-6600; FAX: (615)399-6600

C	rm Apr.08.2008	rm Apr.08.2008
B	rm Feb.19.2008	rm Feb.19.2008
A	rm Nov.20.2007	rm Nov.20.2007
Alter.	Drawn	Checked

Reference	Type	TOPLINE40/B DT-NDV Thermal disinfection, Internal Steam Generation	
	Proj. no.		
	Drawing no.		A Index

Scale: 1:20	Drawn: DB	Date: MAY.13.2011	Checked: MF	Date: MAY.13.2011
-------------	-----------	-------------------	-------------	-------------------